Assessing the performance of an oral biotherapeutic delivery system (OBDS) using intra-duodenal endoscopy delivery in Yucatan minipigs



Shaoying Nikki Lee¹, Cheryl Stork¹, Jeff Smith¹, Bryan Smith¹, Nelson Quintana¹, Chris Wahl¹, Sharat Singh¹ 1. Biora Therapeutics, Inc., CA, USA

INTRODUCTION

Biologics/peptides/nucleic acids are highly effective drugs; however, oral delivery of these therapeutics has proved to be difficult due to the harsh conditions of the upper gastrointestinal tract (GIT) and the poor absorption rate in the small intestinal mucosa. The current state of-the-art technology for a successful oral protein delivery provides around 1% bioavailability when delivered as an oral tablet (Rybelsus® oral Semaglutide).

We aim to develop an oral biotherapeutic delivery system (OBDS) that prevents drug degradation in the upper GIT and increases bioavailability via submucosal injection. The OBDS capsule operates autonomously and provides a needleless injection to deposit the liquid drug payload into the submucosal space of the proximal small intestine.

OBJECTIVE

To develop an intraduodenal endoscopic placement method to place a semi-autonomous OBDS device into the small intestine of swine to allow natural transit, triggering, and submucosal injection for better human translation.

METHODS

In-vivo evaluation of the Pharmacokinetics of PGN-OB1 via Intraduodenal endoscopy placement in the Yucatan Swine

Intraduodenal (ID) endoscopy placement of the OBDS device(s):

• OBDS capsules, which were filled with India ink for the *in-vivo* ink deposition test, or with a variant of adalimumab (PGN-001) for evaluating the pharmacokinetics of PGN-OB1, were attached to the endoscope via the working channel and capsule endoscope delivery device and inserted orally into fasted animals under anesthesia (Figure 2).

- The device was advanced past the pyloric sphincter and expelled from the capsule placement instrument in the proximal small intestine. The capsule then transited naturally and autonomously triggered in the GIT.
- Blood samples post-ID dosing were collected to evaluate the injection efficiency of the OBDS compared with the IV control group (Table 2).



RESULTS

In-vivo OBDS ink deposition in swine

- OBDS ink capsules were placed in proximal small intestine and manually triggered to confirm deployment. (Figure 4A-4C)
- At 24 hr. post-dose, the animal was sacrificed and ink deposition was observed in the small intestine (Figure 4D)



FIGURE 4. Endoscopy placement of OBDS ink capsule In vivo. **A.** placement of OBDS capsule in the duodenum; **B.** manual triggering of OBDS; **C.** Close-up look at the tissue ink deposition; **D.** Ink deposition at terminal necropsy of swine duodenum at 24hr post-deployment.

PRECLINICAL MODEL

- Although the canine is a preferred model for oral therapeutic evalution (as described below), anatomical differences between the canine and human small intestine make it suboptimal for the evaluation of intestinal injection (Table 1).
- A Yucatan minipig model was chosen to better represent the pharmacokinetic properties of submucosal injection in humans. Ex-vivo testing results showed similar tissue ink deposition to human (see poster #105).

TABLE 1. Physiological and Anatomical Difference Comparison

	Gastric pHª	Duodenum pHª	Gastric Emptying Time, Fasted (hr)	Gastric Emptying Time, Fed (hr)	Small Intestinal Transit Time (hr)	Small Intestinal Volume (ml)ª	Small Intestinal Villi Shape
Human	0.4-4 (fasted) 2-4.5 (fed)	5-7	0.66-1 ^b	2-5 ^g	2-4ª	212±110	Fingerª
Swine	1.4-4 (fasted) 4.4 (fed)	6	Variable; 1.4 and up to 20 days ^{c,d}	Variable ^{h,i}	Variable 3-4; 1-2 days ^{d,h}	476±253	Finger ^a
Canine	1.5 (fasted) 3-5 (fed)	6.2	0.4-1 ^{e,f}	Variable; 12-13 ^j	2-3 ^j	300	Long and slender ^a

^aHatton et al. 2015; ^bWorsoe et al. 2011; ^cDavis et al. 2001; ^dHossain et al 1990; ^eMahar et al 2012; ^fKolziek et al 2019; ^gLee et al 2014; ^hGregory 1990; ⁱTreacy 1990; ^jLinbury et al 2012

Preclinical Animal Model Selection for Human Translation

Swine model for proof of mechanism of action and performance test for a tethered device and a semi-autonomous device (see figure 1)

• Similar anatomical and histology features in the small intestine Pros • Better representation of PK for submucosal injection for human translation FIGURE 2. Endoscope Delivery Device for ID placement of OBDS

TABLE 2. Study Design

Group	Test Article	Dose Route	Ν	Nozzles	Dose	PK Timepoints			
1	PGN-OB1	*ID- Uncoated Trigger	13	2 or 6 nozzles	~75 mg	Pre-dose, 0.5, 1, 2, 4, 6, 8, 10, 24, 48, 72, 96, 144, hours post dose (ID)			
2	PGN-OB1	**IV	1	N/A	~75 mg	Pre-dose, 0.25, 0.5, 1, 3, 8, 10, 24, 48, 72, 96, 144, hours post dose (IV)			
*ID: Intraduodenal endoscopy placement; **IV: Intravenous injection									

Pharmacokinetics of PGN-OB1 via endoscopy placement in swine

- All OBDS capsules were successfully advanced through the pyloric sphincter, without early deployment, and were released in the proximal duodenum to naturally transit and deploy in vivo.
- Eight animals showed detectable drug levels (Figure 5), and an oral bioavailability average of 25% (range from 7-55%), excluding an animal showing a late deployment at 72hr post-dose.



FIGURE 5. Plasma concentration of PGN-001 treated with ID and IV over time



 Variable GI transit → higher variability Cons • Prolong gastric emptying time \rightarrow cannot fully evaluate autonomous trigger

Canine model for repeatability and consistency of fully autonomous device

- Similar GI transit and motility to human
- **Pros** Consistent and controllable gastric emptying
 - Ease of oral dosing and repeat dosing for consistency test
- **Cons** May underestimate the bioavailability due to less injection volume/deposition (see poster #105)

Limitations

Swine represents a good model to understand the potential human pharmacokinetics (PK) of submucosal injection, however, variability is expected with an autonomous trigger device due to variable small intestine transit time, motility, gas, and water pockets when compared to the human or canine model (Table 1).

• An OBDS protoype device with tethered triggering in the fixed location in the proximal small intestine showed similar bioavailability but less inter-animal variability ($\sim 26\% \pm 7\%$). Consistency and repeatability of a fully autonomous device will be further examined

FIGURE 1. Preclinical Model Selection

Selected References

• Hatton GB, Yadav V, Basit AW, Merchant HA. Animal Farm: Considerations in Animal Gastrointestinal Physiology and Relevance to Drug Delivery in Humans. J Pharm Sci. 2015;104(9):2747-2776.

• Hossain M, Abramowitz W, Watrous BJ, Szpunar GJ, Ayres JW. Gastrointestinal transit of nondisintegrating, nonerodible oral dosage forms in pigs. Pharm Res. 1990;7(11):1163-1166.

• Koziolek M, Grimm M, Bollmann T, et al. Characterization of the GI transit conditions in Beagle dogs with a telemetric motility capsule. Eur J Pharm Biopharm. 2019;136:221-230.

CONCLUSION

In this study, we have demonstrated that PGN-OB1 can achieve as high as 55% bioavailability of a variant of adalimumab, which is a magnitude greater than current oral protein or peptide delivery technology in the market, and at levels much closer to the subcutaneous route of administration estimated in human trials.

Scan the code to view an electronic version of this poster, or visit bioratherapeutics.com/publications

Poster #106

Presented at Controlled Release

CONTACT INFO

Chris Wahl, MD, MBA VP, Strategy and Business Development Society Annual Meeting, July 13-14, 2022 chris.wahl@bioratherapeutics.com

© 2022 Biora Therapeutics, Inc. All rights reserved. Biora Therapeutics™ is a trademark of Biora Therapeutics, Inc.

bioratherapeutics.com

in the canine model.